

NITRILE GLOVE SIZE CHART

Place your right hand on the sheet making sure that the blue circle fits exactly between your thumb and forefinger like in the image. Keep your fingers together so that your hand spreads out correctly and gently press your hand on the sheet. Measure your hand at the widest part and that is your nitrile glove size.

Disposable nitrile gloves should fit snugly on your hand without feeling like they restrict your range of motion.

A good fit ensures that the glove feels like a second skin.

You should be able to move your fingers normally without over-stretching the glove as you move and bend



Mediqapp BV[©]

SMALL

MEDIUM

LARGE

EXTRA LARGE

<https://Mediqapp.eu>